

**WEAVER**  
and company

565 Nucla Way, Unit B  
Aurora, Colorado 80011 USA

# Nuprep®

Skin Prep Gel



**4oz (114g) Tubes**

Sold as a 3-pack  
ITEM #10-30



**25g Tubes**

Sold as a 6-pack  
ITEM #10-61



**25g Tubes**

Sold as a 20-pack  
ITEM #10-65



**Single Use Cups**

Sold as a 24-pack  
ITEM #10-91

## Target Markets

- EEG Examinations
- Sleep Diagnostics
- Evoked Potential Procedures
- Intraoperative Monitoring
- Audiology
- Biofeedback
- Cardiac Stress Tests
- Holter Monitoring
- Cardiac Rehabilitation Monitoring

## Contact Us

(800) 525-2130

[www.weaverandcompany.com](http://www.weaverandcompany.com)  
[info@weaverandcompany.com](mailto:info@weaverandcompany.com)

## About Nuprep Skin Prep Gel

For quality tracings, proper skin preparation at every electrode site is crucial.

Nuprep Skin Prep Gel gently abrades and moistens the top layer of skin. This lowers the skin impedance with minimal skin irritation and patient discomfort.

Nuprep virtually eliminates the undesirable effects of diaphoresis and unwanted artifact, thus providing the best opportunity to achieve accurate neurodiagnostic and cardiologic tracings.